

June 12, 2002



Sandia's μ ChemLabTM in Industrial Applications

Alex Robinson, Wayne Einfeld,
Steven Showalter, Art Pontau, Bruce Kelley

Acknowledgements:

SNL: Richard Kottenstette, Ron Manginell, Matthew Blain, Pat Lewis, Joy Byrne, Sherri Zmuda, Mathew Moorman, James Sanchez, Fred Gelbard, Dan Horschel, Richard Cernosek

Thermo-ONIX: Jerry Evans, Wallace Trochesette, Jim Bittikoffer

Phillips Petroleum: Joe Cross, John Hays, Chris LaFrancois

Sandia is a multiprogram laboratory operated by Sandia Corporation, a Lockheed Martin Company, For the United States Department of Energy under contract DE-AC04-94AL85000.

S&C FY02





Sandia's Enabling Capabilities Produce Miniature Sensors, Processors, and Communication Systems



Sense, Process, Communicate




**Micromachines,
Microelectronics**



Microsensors



**Photonics,
Microwave Circuits**



**Microelectronics
Development Laboratory
MDL**



**Over 30,000 ft² of
clean room, 0.6 μ m CMOS
Fabrication Facility**

**Integrated Materials
Research Laboratory
IMRL**



**Materials Fabrication and
Characterization, including Plasma
Deposition and Surface and
Interface State Characterization**

**Compound Semiconductor
Research Laboratory
CSRL**



**MOCVD, MBE Deposition,
Electron Beam Lithography,
Reactive Ion Beam Etching**

μ ChemLab Application Strategy


**Sensitive
Selective
Fast**



**Low Power
Hand Held
Low Cost
Versatile**

S&C FY02

Non-proliferation



Counter Terrorism



Military (CW/BW)



Biomedical Diagnostics



Industrial Processes



Environmental



Industrial Hygiene



Food and Water Safety



Microfabricated Components for Sandia's μ ChemLab

Preconcentrator/
Thermal Desorber

4 SAW
Array




One Meter
GC Column

S&C FY02

Broad Applicability

μ Columns Tailored to Specific Applications

- Chemical warfare agents/simulants
- Pharmaceutical drying, VOCs
- Biological analysis, FAMES
- Light gases: CO, CO₂, C₁-C₄



S&C FY02



Ethylene Refinery

Need quantitative measurements of:

- ethane, ethylene, acetylene, carbon monoxide
also present: carbon dioxide, methane



Changes from Sandia gas phase μ ChemLab:

- High, positive pressure
- Diaphragm micro-injection valve
- PDID (pulsed discharge ionization detector)
- permanent installation in hazardous environment




S&C FY02

Current OIT Project

Petroleum Industry: Ethylene Production

Problem Statement:

- Ethylene: high-volume domestic and worldwide production
- 1999 US production ~26 M ton
- Production via batch dehydrogenation of ethane feedstock
- Acetylene is undesirable byproduct





Issues



- Acetylene content controlled by catalyzed reaction in ARU requiring rapid process feedback control
- Disrupt conditions in ARU occur on short ($< 3\text{min}$) time scales
- Reactor runaway requires dumping to flare (\$40K/hour)
- Current process methods are expensive to implement and maintain

Goals

- Reduce product loss
- Increase yield efficiency
- Reduce overall energy consumption
- Lower instrumentation costs



Advantages of Integrated Microsystem




- Replace costly instruments (\$30K/each) in expensive shelters (\$300K/each) with low-cost, MEMS technology
- Small-footprint, on-line process monitoring (intrinsically safe)
- Faster (1 min cycle time) process monitoring
- Improve real-time process control for enhanced process energy efficiency
- Allow increased sensor density for ARU stages

S&C FY02

Estimated Savings

- Faster cycle-time analyzers will reduce product loss from 12.6 M lbs. to 2.5 M lbs.
- Equivalent to \$400K/year/ARU product cost + energy efficiency savings

Present Ethylene Product Yield *



*Data from a Phillips ARU over a 16-mo. interval



Budget

DOE-OIT	FY 2001	FY 2002	FY2003
Sandia National Labs	\$300K	\$325K	\$325K

Industry in-kind

Thermo-ONIX	\$275K	\$275K	\$225K
Phillips Petroleum	\$25K	\$50K	\$100K

Total



	\$600K	\$650K	\$650K
--	--------	--------	--------

3 Year Total

\$1900K

Technical Accomplishments

- Demonstrated success of micro-packed column to separate gases of interest
 - 20 cm long channel; 300 x 300 μm cross section; packed with Supelco carboxen 1000E, 75-100 μm
- Identified micro-PDID detector as best candidate for system
 - evaluation in progress
- Completed micro-valve design
 - preliminary fabrication in progress




S&C FY02




On-Chip, Packable Micro Column

- Enables use of existing column packing materials and technology
- Allows necessary retention of light molecules for adequate separations
- Improves resolution
- Patent filed



S&C FY02


Micropacked Carbon Column Effectively Separates Acetylene in Ethylene



10 μ L injection; 60 °C isothermal; 10 p.s.i.g. N₂ Supelco carboxen 1000E packing

S&C FY02

Valco micro Pulsed Discharge Ionization Detector



commercial diaphragm valve/ μ GC/PDID


Advantages:

- Allows CO, CO₂ detection
- Flame-proofing not required
 - instrument located at-line
- Sensitive
- Robust



S&C FY02

Microfabricated Diaphragm Sample Valve Design (in experimental stages)



Advantages:

- Lower production cost
- Integrates system functionality
 - smaller dead volumes
- Robust (?)

- Pyrex layer for ports
- Silicon layer with analyte conduits and valve features
- Kapton or silicon diaphragm
- Silicon layer with P-actuation conduits

S&C FY02



Important Project Milestones

Over Project Lifetime

- Identify best column configuration - completed
- Identify best micro-detector – completed
- Design micro-gas sampling valve - completed
- Fab micro-gas sampling valve – in process
- Optimize micro-packed column – in process
- Optimize micro-PDID detector – in process
- Calibration methodology - pending
- Develop system architecture - completed
- Package components – pending
- Field test – pending
- Explore commercialization pathways – in process


S&C FY02

Tailored Instrument Delivered to a Customer



Temperature Ramp with TEC





98 cm x 75 μm x 400 μm open μGC



Applications in fermentation, synthesis, and product drying in pharmaceutical production.

S&C FY02

Sandia Liquid Phase μ ChemLab™



S&C FY02



Global Conclusions

- **μChemLab systems can meet many important applications.**
- **Microfabrication provides small size, low power, rapid heating, low dead volumes, and potentially low cost.**
- **We have designed, fabricated, and demonstrated:**
 - **Rapid, low power sample preconcentrators; use of injection valves**
 - **Miniature gas chromatograph columns (open and packed)**
 - **High sensitivity sensors (SAW, PDID, TCD)**
- **Small, autonomous chemical analyzers are currently being assembled and tested for various applications**
 - **government:** military, 3-letter agencies
 - **industrial:** petroleum, pharmaceutical, energy production, others?

Short list of Sandia Sensor Technologies


- μ ChemLab




- Micro Ion Mobility Spectrometer




- Chemiresistor and SAW Arrays for pattern recognition




- Fringe Field Sensors




- Ionizing Radiation Sensors



- Fiber Optic Chemical Sensors



- Programmable Diffraction Grating



- Combustible Gas Detector
- Hydrogen Sensors
- Portable Cloud Point Detector

S&C FY02